



cicCartuja

centro de  
investigaciones científicas  
isla de la cartuja

Consejo Superior de Investigaciones Científicas  
Junta de Andalucía  
Universidad de Sevilla

## Specific Use Rules for DLS Service

- ✓ Only authorised users will have access to the equipment. They will be given a user name, a password to access the computer controlling the DLS and a secure folder to save their experiment data.
- ✓ To use the equipment, users must reserve on the Biomolecular Interaction Platform web application in advance. Reservations are personal.
- ✓ During the operation of the equipment, do not change the environmental conditions (light, temperature...). Affects the analysis result.
- ✓ Turn on the equipment half an hour before starting the analysis. The minimum reserve time is one hour.
- ✓ To ensure that there are no delays for the next user, make sure everything in the equipment is back in the correct conditions (e.g. standard temperature 25 °C) when the shift is finished.
- ✓ At the end of the shift, users must logout leaving and the used material in its right place.
- ✓ Users must communicate any incident or observation that is relevant to the person responsible through the email [bip@ciccartuja.es](mailto:bip@ciccartuja.es).
- ✓ Authorised users are responsible for their own consumables.



cicCartuja

centro de  
investigaciones científicas  
isla de la cartuja

Consejo Superior de Investigaciones Científicas  
Junta de Andalucía  
Universidad de Sevilla

- ✓ Users will be responsible for refund if any failure/breakdown is caused by negligence or misuse of equipment or incorrect setup of experiments.