



cicCartuja

centro de  
investigaciones científicas  
isla de la cartuja

Consejo Superior de Investigaciones Científicas  
Junta de Andalucía  
Universidad de Sevilla

## Instructions to gain access to MALDI-TOF Mass Spectrometer

- ✓ Read the *General Rules for Access and Use* at the *General Guidelines* section. You need have the *BIP Project Form*, of your Principal Investigator, and the *BIP User Form* filled, to gain access to the Biomolecular Interaction Platform (find there the links of both fillable PDF forms).
- ✓ The first time you contact us, please send e-mail to the responsible at [bip\(at\)ciccartuja\(dot\)es](mailto:bip@iccartuja.es).
- ✓ Due to the technical qualification required, the equipment will be accessed through the responsible.
- ✓ Fill the *MALDI-TOF Sample Submission Form* with information of samples. Send the completed form and samples to the responsible. Please
- ✓ The amount of material required varies. In general, 10 $\mu$ L of a 10  $\mu$ M solution of sample is sufficient for most analyses. In the case of large molecules, this concentration should be considered a minimum.