



cicCartuja

centro de
investigaciones científicas
isla de la cartuja

Consejo Superior de Investigaciones Científicas
Junta de Andalucía
Universidad de Sevilla

Instructions to gain access to MALDI-TOF Mass Spectrometer

- ✓ Read the *General Rules for Access and Use* at the *General Guidelines* section. You need have the *BIP Project Form*, of your Principal Investigator, and the *BIP User Form* filled, to gain access to the Biomolecular Interaction Platform (find there the links of both fillable PDF forms).
- ✓ The first time you contact us, please send e-mail to the responsible: *Javier López Prados* ([rmn\(at\)cccartuja\(dot\)es](mailto:rmn@cccartuja.es)).
- ✓ Due to the technical qualification required, the equipment will be accessed through the responsible.
- ✓ Fill the *NMR Sample Submission Form* with information of samples. Send the completed form and samples to the responsible. Please
- ✓ The amount of material required varies. In general, 10 μ L of a 10 μ M solution of sample is sufficient for most analyses. In the case of large molecules, this concentration should be considered a minimum.