

Specific Use Rules for MALDI-TOF Mass Service

- ✓ Due to the technical qualification required, the equipment will be accessed through the responsible. Please contact directly to the person in charge or by e-mail to bip@ciccartuja.es.
- ✓ Fill the MALDI-TOF Sample Form with information of samples. Send the completed form and samples to the responsible.
- ✓ The amount of material needed varies. In general, 10 μ L of a 10 μ M solution of sample is sufficient for most analyses. In the case of large, complex molecules, the concentration should be considered as a minimum, and ideally close to 100 μ M.